



SLIM^{ex}

NucleOSS Dental Implant System **SURGICAL GUIDE**



This manual contains recommended instructions for use for surgical applications of the NucleOSS Dental Implant system. Treatment and surgical planning should be based on clinical data and the physician's evaluation of the case.

FEBRUARY 2024 EN

Dear Clinician

Şanlılar Ltd.Şti., the manufacturer of the NucleOSS brand, was established in Izmir in 2001 and produced Turkey's first dental implant. In the light of meticulous R&D studies and knowledge, we continue to produce more innovative and user-friendly products for you with the excitement of the first day, by adopting world-class sensitive production protocols. In order to increase the quality of the operation and to enable more comfortable treatment processes, we design and develop quality and innovative products, and offer them to you.

We established the TFI (Together for Implantology) academy in order to create a platform where dentists can share their ideas, studies and information about implants within an ethical framework and thus contribute to the development of the field of implantology. As NucleOSS, we adopt reliable, quality and customer-oriented production with transparent, ethical and principled management.

We work more and more every day within the framework of our open to change and innovative approach, which is the first condition of progress and growth. We design new products, constantly change and develop with our understanding of management and production that values people.

We are growing thanks to your support... We are enchanting...



Salih ŞANLI


PRESIDENT

| <u>Topics</u> | <u>Page No.</u> |
|--|-----------------|
| SlimEX Implant - Bone Level Implant | 3 |
| Implant Features | 4 |
| Diameter and Length Options | 5 |
| Package Content | 6 |
| Cement-retained Solutions | 7 |
| Custom Solutions | 8 |
| Screw-retained Solutions | 9 |
| Overdenture Solutions | 10 |
| Preoperative Planning | 11 |
| Pre-Operational Procedures | 12 |
| T6 Standard and T6 TorQ Surgical Kits | 13 |
| Drill Protocol | 14 |
| Drill Protocol - Auxiliary Parts | 15 |
| Opening the Implant Package | 16 |
| Carrying the Implant | 17 |
| Implant Placement | 18 |
| Technical Specifications of T6 Surgical Kits | 19 |
| Initials Drills | 20 |
| Pilot Parallel Pins | 21 |
| Hex Drivers | 22 |
| Ball(O-ring) and Screw-Retained Tightening Parts | 23 |
| Depth Gauge - Screw Driver - Ratchet | 24 |
| Tightening - Torque Values | 25 |
| NudeOSS Dental Implant System - Implant Prospectus | 27 |
| Notes | 30 |

SLIMex

Offers thin and
durable solutions
in the crests where
bone thickness is
insufficient.

Bone Level Implant



extra slim **extra durable**

Different and wide
prosthetic feature, gingiva
height abutment option.



Slimex implant the newest member of the NucleOSS product family, offers thin and durable solutions in the crests where bone thickness is insufficient.

The certified unique Maxicell surface provides a safe and high survival rate for Slimex implants.

Slimex, which has a durable and reliable structure thanks to its external design, satisfies all implant treatments with its unique and wide superstructure options.

NucleOSS offers the Slimex implant, developed with 20 years of R&D experience and the support of specialist physicians, to the use of you, precious physicians.

► Cemented Solutions



► Overdenture Solutions



► Screw-Retained Solutions




► Digital CAD/CAM Solutions



► Custom Solutions



> Hexagonal External Connection



Thanks to hexagonal external connection of NucleOSS Slimex Implant, abutments can be applied with 6 different positions within 60° rotation.

> Double Thread

Double Thread and reverse buttress form provide faster and easy loading. Increases contact area between bone and implant.

> Strong Side Walls

> Ø 3.0 Diameter


3.0 mm implant diameter for narrow crests.

> Cutting Edges

It enables the implant to move into the implant bed steadily and without damage when placing the implant.


> Maxicell Surface

Proven and certificated clean and hydrophilic unique surface.




Certified Maxicell Surface


Zoom
x100




Zoom
x1000



Zoom
x1500



Zoom
x3000







Slim^{ex} bone level narrow implant designed with modern and aesthetic details; It offers a comfortable surgical process to the physician and the patient with its safe surgical protocol. It brings bold smiles with its wide superstructure range. It refreshes your self-confidence with its smiling designs and advanced strength.

| Slim^{ex} implant is provided in four different sizes.


Ø3.0




| Lenght | Product No |
|----------|------------|
| L : 08 > | TS 3008 |



| Lenght | Product No |
|----------|------------|
| L : 10 > | TS 3010 |



| Lenght | Product No |
|----------|------------|
| L : 12 > | TS 3012 |




| Lenght | Product No |
|----------|------------|
| L : 14 > | TS 3014 |

Package Content

SLIM^{ex}

Slimex implant is provided in blister and sealed sterile tube. The implant package includes; comfort cap, impression post, abutment screw and standard abutment GH:2mm.




Slimex implant boxes are opened by tearing over the line that created for security purposes.

The first safety barrier of Slimex implants are thermoforms. The package is opened by holding the opening end on the edge and the tube is taken out.

Package content consists of 3 parts.




(extra)




| Product Image | Material | Prosthetic Platform Ø (mm) | Gingival Height GH (mm) | Description | Code |
|---------------|----------|----------------------------|-------------------------|-------------|------|
|---------------|----------|----------------------------|-------------------------|-------------|------|


Standard Straight Abutments

| | | | | | |
|---|--------|-------|-----|----------|----------|
|  | Ti Gr5 | ● 4.2 | 0.6 | Straight | TS 42003 |
| | | | 1.6 | Straight | TS 42004 |
| | | | 2.6 | Straight | TS 42005 |
| | | | 3.6 | Straight | TS 42006 |

Angled Abutments

| | | | | | |
|---|--------|-------|---|-----|----------|
|  | Ti Gr5 | ● 4.2 | 1 | 15° | TS 42007 |
| | | | 2 | 15° | TS 42008 |
| | | | 3 | 15° | TS 42009 |
| | | | 4 | 15° | TS 42010 |

Auxiliary Parts

| | | | | | |
|---|-----------------|-------|---|-------------------------------|----------|
|  | Plastic | ● 4.3 | | Impression Cap | TO 32909 |
| | Ti Gr5 | | | Abutment Screw | TS 18000 |
| | | | | Open Tray Impression Screw | TS 18002 |
| | Ti Gr5 | ● 4.2 | 3 | Open Tray Impression Post | TS 32606 |
| | | | 3 | Close Tray Impression Post | TS 32603 |
| | Plastic | | | Impression Cap for Close Tray | TO 32912 |
| | Stainless Steel | ● 3.5 | | Analogue | TS 32206 |
| | | | | | |

The abutments are made of TiGr5 material.

The abutments are provided with a connecting screw.

The basis suitable for planning can be selected from the table.



Pack Content




Pack Content


Cemante offers solutions in different prosthetic platform diameters and gingival heights for prosthetic solutions. Angled abutments offer solutions in different prosthetic platform diameters and gingival heights for cemented angled prosthesis solutions. It has 15° angle options.

| Product Image | Material | Prosthetic Platform Ø (mm) | Gingival Height GH (mm) | Description | Code |
|---------------|----------|----------------------------|-------------------------|-------------|------|
|---------------|----------|----------------------------|-------------------------|-------------|------|


Casting Caps

| | | | | | |
|---|--------|-------|--|-----------------------|----------|
|  | Delrin | ● 3.5 | | Burn Out Cap (NO HEX) | TS 35000 |
| | | | | Burn Out Cap (HEX) | TS 35001 |
| | Ti Gr5 | | | Abutment Screw | TS 18000 |



Auxiliary Parts

| | | | | | |
|---|-----------------|-------|--|----------------------------|----------|
|  | Ti Gr5 | ● 4.2 | | Open Tray Impression Screw | TS 18002 |
| | | | | Open Tray Impression Post | TS 32606 |
| | | | | Close Tray Impression Post | TS 32603 |
| | Plastic | | | Impression Cap | T0 32912 |
| | Stainless Steel | ● 3.5 | | Analog | TS 32206 |

CAD-CAM Abutments

| | | | | | |
|---|--------|-------|-----|--------|----------|
|  | Ti Gr5 | ● 4.2 | 0.6 | NO HEX | TS 42012 |
| | | | 2.4 | NO HEX | TS 42013 |
| | | | 2.4 | HEX | TS 42014 |
| | | | 0.6 | HEX | TS 42015 |








CAD-CAM Auxiliary Parts

| | | | | | |
|---|-----------------|-------|--|------------------------------|----------|
|  | Ti Gr5 | ● 4.8 | | T6 Multi Unit Digital Analog | T6 32210 |
| | Stainless Steel | ● 3.5 | | Slimex Digital Analog | TS 32209 |
| | Ti Gr5 | | | Digital Analog Screw | T0 2015 |
|  | PEEK | ● 3.5 | | Scan Post | TS 42011 |
| | Ti Gr5 | | | Abutment Screw | TS 18000 |



CAD/CAM Ti-base abutments are designed for physicians and laboratories making dental designs on digital platforms. Hex connection options are available for single tooth treatments and No Hex connection options are available for multiple tooth treatments.

Universal Casting Caps are designed for physicians and laboratories who design personalized superstructures with the casting technique. It is made of Delrin material. Hex (White) for single tooth and No-Hex (Black) for multi-tooth applications are available. Supplied with one cast cap and abutment screw.

| Product Image | Material | Prosthetic Platform Ø (mm) | Gingival Height GH (mm) | Description | Code |
|---|----------|---|----------------------------|-------------------------|------------------------|
| Multi Unit Abutments (0°) | | | | | |
|  | Ti Gr5 |  4.8 | 1.5 | Straight | TS 47000 |
| | | | 2.5 | Straight | TS 47001 |
| | | | 4.5 | Straight | TS 47002 |
|  | Ti Gr5 |  4.8 | 4.4 | Comfort Cap | T6 32915 |
| | | | 5.7 | Comfort Cap | T6 32921 |
|  | Ti Gr5 |  4.8 | | Occlusal Screw | T6 16009 |
| Multiple Teeth (NO HEX) | | | | | |
|  | Ti Gr5 |  4.8 | | Impression Cap | T0 32912 |
| | | | | Impression Post (close) | T6 32614 |
| | | | | Impression Post (open) | T6 32613 |
|  | Ti Gr5 |  4.8 | | Ti-Base Abutment | T6 32813 |
| | | | | Temporary Titanium Cap | T6 32916 |
|  | PEEK |  4.8 | | Scan Post | T6 32037 |
| | Delrin | | | Burn Out Cap | T6 32917 |
| | Ti Gr5 | | | Analog | T6 32207 |
| Single Tooth (HEX) | | | | | |
|  | Plastik |  4.8 | | Impression Cap | T0 32912 |
| | Ti Gr5 | | | Impression Post (close) | T6 32617 |
| | | | | | Impression Post (open) |
|  | Ti Gr5 |  4.8 | | Ti-Base Abutment | T6 32814 |
| | | | | Temporary Titanium Cap | T6 32919 |
|  | PEEK |  4.8 | | Scan Post | T6 32038 |
| | Delrin | | | Burn Out Cap | T6 32918 |
| | Ti Gr5 | | | Analog | T6 32208 |



Pack Content





Pack Content



Multi unit abutments are shipped with NO HEX option, green color burn out cap and two occlusal screws. Multi unit abutments are shipped with HEX option, blue color burn out cap and two occlusal screws.

| Product Image | Material | Prosthetic Platform Ø (mm) | Gingival Height GH (mm) | Description | Code |
|---------------|----------|----------------------------|-------------------------|-------------|------|
|---------------|----------|----------------------------|-------------------------|-------------|------|


Ball Abutment

| | | | | | |
|---|--------|---|-----|--|----------|
|  | Ti Gr5 |  4.5 | 1.5 | | TS 45003 |
| | | | 2.5 | | TS 45004 |
| | | | 4.5 | | TS 45005 |




Ball Abutments Auxiliary Parts

| | | | | | |
|---|-----------------|---|--|-----------------------------------|------------|
|  | Stainless Steel |  3.5 | | Analog | TO 4521 |
| | Ti Gr5 | | | Impression Cap | TO 32615 |
| | | | | Metal Ring | TO 5200 |
| | Silicon | | | Silicon Ring | TO 4500-3 |
| | | | | Silicon Ring (for impression cap) | TO 4500-18 |

Equator® Abutments



| | | | | | |
|---|--------|--|-----|--|----------|
|  | Ti Gr5 | | 1.5 | | TS 45000 |
| | | | 2.5 | | TS 45001 |
| | | | 4.5 | | TS 45002 |

Equator® Abutments Auxiliary Parts

| | | | | | |
|---|---------|---|--|---------------------------------|----------|
|  | Ti Gr5 |  4.5 | | Impression Post | 044 CAIN |
| | | | | Analog | 144 AE |
| | | | | Tightening Part (Ratchet) | 774 CHE |
| | | | | Tightening Part (motor) | 760 CE |
| | | | | Cap Inserter and Extractor Tool | 485 IC |
|  | Ti Gr5 | | | Metal Housing | CA008 |
| | | | | Transparent/Standart | CN051 |
| | | | | Pint/Soft | CN052 |
| | Silicon | | | Black/Laboratory | CN068 |

**Pack Content****Pack Content**

Ball abutment solutions are offered to users with different gingival heights in removable total prosthesis treatments. It is made of TiGr5 material. Equator® abutment solutions are offered to users with different gingival heights in removable total prosthesis treatments. It is made of TiGr5 material.



While planning the superstructure in dental implant - prosthesis treatment, the positioning of the implant in the patient's mouth constitutes the basic structure.




Implant abutments should always be planned axially and coaxiality should be provided with the lateral teeth. The anatomical condition of the patient should be considered when deciding on the diameter, length, model and number of implants and abutments.

Below are the recommended minimum placement sizes when planning the Slimex implant in the patient's mouth.

Diameter & Length Selection Suggestions :

The image shows the recommended areas of use according to the diameter of the Slimex bone level implant. For treatment and implant selection, the clinician should evaluate the patient, the case and the intended restoration. Use outside the specified regions is the choice and responsibility of the clinician.

Implant Positioning:



In order for the implant to provide primary stability in the bone, there should be a minimum of 1.0 mm bone thickness on the palatal and vestibule surfaces.

In cases where the bone thickness is below 1.0 mm, the augmentation method can be applied.

Pre-Operational Examination

General evaluation of the patient should be made with clinical and oral radiographic tests, and mucous membrane, jaw morphology, dental history, prosthetic history and oral dysfunctional signs should be carefully examined before the operation. Radiographic analyzes should be used in the examination of bone topography.

In addition to clinical examination, preliminary radiographic evaluations are a basis for determining the patient's suitability for implant treatment.

If the patient is found suitable for treatment, a more detailed clinical examination should be performed for the treatment area and the countervailing area. Any pathological condition detected in the jaw should be treated before the implant operation.

Pre-Operation Planning

While performing the planning before the operation, the final prosthetic treatment should be considered and planning should be done accordingly. First of all, the intended final prosthetic treatment should be planned, and the type and number of implants should be decided accordingly. The prosthetic material to be used should be chosen according to the region where the implant will be placed. While planning the treatment, all stages from the healing process to the final restoration should be evaluated with the intended final treatment in mind.

Articulator-mounted models provide information about the connection between the tooth and the jaw. Wax modeling from the edentulous area makes an important contribution to correct planning. The most appropriate treatment planning can be achieved by examining the occlusal picture and evaluating factors such as load distribution and preferred slots for implants.

The Transparent Radiographic Implant Guide, showing the magnified dimensions of the implants at different rates, will contribute to the most appropriate selection of the point, direction and length of implant.

During the surgical operation, factors such as the primary stability of the implant, single-stage or two-stage surgery, immediate implant placement after extraction, and expected recovery time before loading should be considered as well as the intended final treatment approach.

The recommended healing period before loading is 12 weeks, unless there is a different medical decision taken by the physician.


In each case, the loading time for each implant should be determined by the clinician by carefully examining bone quality, quantity, density, primary stability, restoration style, and loading conditions.

Before starting the treatment, hole process should be clearly explained to the patient by specifying the results of the examination performed before the operation, the required treatment method, the expected result after the treatment, the required care and possible risks.

Bone Classification


Since the thickness and density of the bones of every patient vary, implant placement and surgery also varies for each case. According to the case data identified at the examination, the appropriate route should be followed for loading and implant surgery. Therefore, the NucleOSS surgical set and protocol were prepared taking into account these functional needs.

D1




Dense cortical
Anterior Mandible

D2



Porous cortical and thick trabecular
Anterior Mandible Posterior
Mandible Anterior Maxilla

D3



Porous cortical and (thin)
tight trabecular
Anterior Maxilla
Posterior Mandible
Posterior Maxilla

D4



Tight trabecular
Posterior Maxilla

The NucleOSS Dental Implant System has been developed considering the natural oral anatomy and the intended final prosthetic restoration.

Positions of implants are recommended based on crown-implant compatibility. However, while deciding on the placement of the implant, the size and shape of the implant, soft tissue healing, support required for the intended final restoration, bone volume, and connection to the neighbor tooth should be considered.

Based on mechanical considerations, it is always recommended to use the largest diameter implant possible. This recommendation should be especially taken into account for the rear areas where the loading forces are higher. In all cases, loading conditions and the intended final prosthetic restoration should be considered when planning implant surgery and determining the number of implants.

T6 Standard and T6 TorQ Surgical Kits


SLIMex

In order to place Slimex implants in the jawbone, the point and pilot drills in the T6 standard and T6 torq surgical sets should be used.


CRH-109
T6 Standard Surgical Set



CRH-115
T6 TorQ Surgical Set



| | | | | | |
|------------|---|----------|----------|----------|----------|
| Lenght | > | L:08 | L:10 | L:12 | L:14 |
| Product No | > | TO 25043 | TO 25044 | TO 25045 | TO 25046 |



Step 1 : Point drill is done to determine the slot location.

The final drills for SlimEX implants are the standard pilot drills found in the T6 standard and T6 TorQ surgical sets. The drills are with stoppers and there are separate drills for each implant size.


The slot preparation is completed by entering the slot opened with a point drill with the drill selected according to the implant size.

Drills should be run with coolant at max. 800 rpm.


Step 2 : The slot preparation is completed by entering the slot opened with a point bur with the bur selected according to the implant size.

Step 3 : The implant (with the carrying part) is carefully placed in the opened slot.

Note : Drill length is 1 mm longer than implant length. This should be taken into account in surgical planning.



Step 4 : Finally, the necessary healing caps are placed on the implant and left to the osseointegration process.



Attention ! : NucleOSS standard pilot drills should be used in the drill protocol of the implants.

The Drills are designed 1 mm longer than the implant length. Therefore, if the drill stopper comes into contact with the bone, the bur will open the implant socket 1 mm deeper than the implant length.

Drill lengths should be taken into account when planning the drill for the implants.



TO 25024
Point Drill



TO 0103
Round Drill



TO 2012-1
Frez Extension



TO 25320
Pilot
Parallel Pin



TO 25212
Ratchet
Adaptor



TO 2024
Hex Driver
Short (Ratchet)



TO 2025
Hex Driver
Long (Ratchet)



TO 2009
Hex Driver
(Motor)



TO 25251
Hex Driver
(Handheld)



TO 25205
Ball Abutment
Tightening Parts



TO 25206
Screw-Retained
Abutment
Tightening Parts




TO 2005-1
Tightening Parts
Long (Ratchet)



TO 2006-1
Tightening Parts
Short (Ratchet)



TO 2008
Tightening Parts
(Motor)



TO 25312
Depth Gauge



TO 25900
Ratchet



TO 2022
Screwdriver



Opening the package

Slim^{ex} Implant boxes are opened by tearing over the traced line created for security purposes.

Opening the thermoform

The thermoforms, which are the first safety barrier of Slim^{ex} implants, are opened by holding the opening tip on the edge and the tube is taken out.

Opening the tube cap

The tube cover, which is another safety barrier, is opened by turning it counterclockwise. There is a security ring attached to the covers with thin nails. When the cap is turned, the tabs break off and the safety ring is separated from the cap and remains on the tube.

If this ring is separated from the cap at first use, it indicates that the tube has been opened before.


Caution : Thermoform or tube cap unlocked implants should not be used.

The tube is removed by turning it upside down from the closest distance on the carrier glove or sterile clean cloth in the Slim^{ex} implant with the opened cover.




The carrier cover is held with the help of the thumb and forefinger, turning it slowly with the help of the thumb and index fingers of the other hand and pulled upwards.

Caution : Caution: Once the carrier is capped, the implant is released, so the carrier should never be tilted more than 60°. Otherwise, the implant may come off and fall.



After the implant carrier cover is opened, the implant is transported and placed in the area where it will be placed with the help of Slimex carrying parts in T6 surgical kits.



Ratchet


Screwdriver

Motor



Attention : The red silicone ring on the ends of the transport clamping parts prevents the implant from falling during transport. Therefore, make sure that the silicone ring is in place and undamaged before use.






The Slimex implant, which is removed from the carrier, is placed in the patient's mouth as soon as possible without touching anything, into the implant socket opened according to the protocol.

While placing the Slimex implant with a motor or by hand, turn it clockwise at a maximum speed of 15 rpm until the shiny platform of the implant with a diameter of $\varnothing 3.5\text{mm}$ comes into contact with the bone surface.

Caution : Implant placement torque is recommended as 35-40 Ncm, in order to minimize bone damage at torques higher than these values, the implant should be removed and the socket checked.



Note : The O-ring abutment clamping piece is compatible with the external hex connection of the implant and can be used after the implant is placed.



Placement of the comfort cap

Step 1 : After making sure that the implant is placed, the carrying clamping part on the implant is removed by turning it counterclockwise with the help of a manual hex wrench.

Step 2 : The healing head in the second eye of the thermoform is moved with the help of the hex key. Before inserting into the implant, the inside of the implant is cleaned with the help of a probe, then the cap is placed on the implant by turning it clockwise and tightened with a maximum torque of 10 Ncm.




Attention: The healing head supplied with the standard package is 3 mm high. Healing head with a height of 6 mm can be selected for higher gingiva.


Technical Specifications of T6 Surgical Kit

SLIMex

Dimensions



Opening Positions



NudeOSS surgical sets are designed in direct collaboration with healthcare practitioners, giving you the best possible match between the ergonomic design of the surgical tray, the surgical instruments it contains and its size.

The surgical set cover can be removed with one easy move.

Lightweight, easy to use, stackable, compact and roomy. It is easy to carry and guarantees the safety of the tools it contains, prevents them from being misplaced during transport.

The surgical set is made of highly impact resistant plastic materials suitable for autoclave sterilization. Plastic raw material can withstand more than 1000 cycles of steam sterilization.

T6 surgical sets consists of two parts, a movable middle shelf where surgical burs and tightening parts are located, and a lower section where torque ratchet and other hand tools are located.

Thanks to its compact and convenient design, it can be sterilized as one piece without having to be taken the parts inside out of the box in clinical autoclave devices.

Surgical sets are suitable for autoclave sterilization in the B134°C Prion (134°C/2.15 Bar/18 min) cycle.

Set cover can be used in two different positions by being opened 90° or 210°. When the cover is closed, the parts inside are locked automatically.

Surgical kit consists of drills, parallel pins, neck expanders, hex wrenches and hand tools.

Thanks to the flat and embedded type of silicones holding the surgical parts, they are dirt-proof and easy to clean.

Surgical drills are made of high-quality stainless steel, and after heat treatment, they have a long-lasting cutting edge with sharpening techniques.

T6 Standard Surgery Kit



T6
standard
CRH-109

T6 Torq Surgery Kit



T6
torq
CRH-115


General Info




Round Drill - T0 0103



Point Drill - T0 25024



Drill Extender - T0 2012-1



The parts in the T6 surgical sets are used with the micromotor or the NudeOSS ratchet or screwdriver in the surgical kits.

In addition, parts such as connecting screws placed on the implant, implant caps, healing caps or connecting screws placed on the abutment can also be used with the manual hex wrench included in the surgical kits.

Caution!

Drills must be correctly attached to the micromotor and should not be operated until they are properly seated.


Initial drills are used to prepare the cavity during osteotomy.

Round drill is used to make necessary corrections by filling at the beginning of the osteotomy if the bone surface in the area where the implant will be placed is damaged or irregular.


The point drill is used to determine the place where the implant will be placed on the bone and to allow the next drill to advance comfortably.

After the necessary corrections are made with the round drill, the milling is first started with the point drill according to the surgical protocol.

The drill extender is used to extend the drill lengths by approximately 16 mm in cases where the drill length is insufficient.




The distance is increased with drill extender in cases where the micromotor is not able to advance in the mesial or distal of the implant area.




Parallel pins are used to control the parallelism with the adjacent teeth in single implant treatments, and to control the parallelism between the implant cavities during surgical intervention in multiple implant treatments

T6 surgical sets have four parallel pins.

The lower part of the parallel pin is $\varnothing 2$ mm and can be used after point drilling.



Hex Driver - Transfer



Hex wrench is used for fixing healing components, prosthetic and impression parts.

It is made of stainless-steel material.

It can be procured in a surgical set and in one piece.


The tip of the hex wrenches used manually and with micromotor is in hex structure and conical shape.

It can be used for carrying connecting screws and healing caps and for light tightening.

Torque up to a maximum of 10 Ncm can be made with these hex driver.

Note: Final torque must never be done. Because the hex slot of the screw or the tip of the hex wrench may be damaged.

Hex Driver (torque)



The screw is held with a hex wrench and placed on the part it will be fixed to. The hex driver is gently turn clockwise to tighten the screw.

Tightening drivers are used for fixing healing components, prosthesis and impression parts.

It is used with a ratchet

It is made of stainless-steel material.


It can be procured in a surgical set and in one piece

The tip of the hex wrenches used for final torque is straight. It is only used for torque with NucleOSS ratchet. It has no retention, it is not used for carrying.



The screw is tightened by turning the ratchet clockwise. By following the torque values specified at the end point, the desired torque value is set and the tightening process is completed.

O-Ring Abutment Tightening Part (ball)




O-Ring abutment tightening part is used to place ball abutments on a T6 implant.

They are made of stainless steel.

It can be used with a ratchet or screwdriver.


The ball abutment is carefully inserted into the body. The ball abutment tightening part (T0 25205) is used to fix the abutment to the body.



Ball abutment should be torqued to a maximum of 30 Ncm using the NucleOSS ratchet.



Screw-Retained Tightening Part



Threaded abutment tightening part is used to place multi-unit straight abutments on T6 implants.

It is made of stainless steel.

It can be used with a ratchet or screwdriver.


Multi-unit straight abutment has a threaded connection structure. When placing the abutment on the implant, attention should be paid to the tightening and torque values.

Caution!


Multi-unit straight abutment should be torqued to a maximum of 30 Ncm when being placed on the abutment!

If the torque value is less or more than the specified value, loosening or breakage of the abutment may occur.

Multi-unit angled abutment is carefully placed on the implant by holding the carrying part and fixed to the implant with the Connecting screw included in the package. The connecting screw should be torqued to a maximum of 30 Ncm with the help of a NucleOSS torque ratchet.



Depth Gauge - TO 25312




It is one of the items that are included in both the T6 Standard and T6 TorQ surgical kits.



In order to be sure of the depth of the opened slot(cavity), a control should be performed using a depth gauge.

Screwdriver - TO 2022



It is one of the items that are included in both the T6 Standard and T6 TorQ surgical sets.

It is offered as an alternative to the NucleOSS ratchet for the use of implant carrying/tightening parts, Hex wrenches, threaded abutment tightening part, ball and equator tightening parts.



Torque Ratchet - TO 25900






The torque value of the NucleOSS ratchet is from 0 Ncm to 70 Ncm.

When the needle is at the top, it shows about 80 Ncm.

The ratchet is used by NucleOSS to tighten all the abutments and screws at the torque values specified in the product catalogs.

It is made of stainless-steel material.

It can be procured in a surgical set and in one piece.



Using as Torqueless Ratchet;

Hold the torque arm and torque needle together. Press down on the wheel with your finger to hold the ratchet in place.

Caution! Direction arrow should point to the top during tightening and to the bottom during untightening.

Using as Torque Ratchet:

Press down on the wheel with your finger to hold the ratchet in place. Then pull the torque arm lightly in the direction of the arrow until you reach the desired torque value.










Caution! The torque arm should not be turned more than the torque values specified on the ratchet. Otherwise, the torque will be incorrect. Direction arrow should point to the top during tightening and to the bottom during untightening

Using Ratchet for Removal:

Use the ratchet with the direction arrow facing down.



| Products | Tighten Parts | Torque Values |
|--|--|---|
| <p>Comfort Caps</p>  |  <p>T0 25251 Manual Hex Driver</p> | <p>Manuel/ Ratchet (max. 10 Ncm)</p> |
| <p>Cemented Abutment Screws</p>  |  <p>T6 2024 Short (ratchet)</p> <p>T6 2025 Long (ratchet)</p> | <p>Max. 30 Ncm</p> |
| <p>Cad/Cam Ti-Base Abutment Screws</p>  |  <p>T6 2024 Short (ratchet)</p> <p>T6 2025 Long (ratchet)</p> | <p>Max. 30 Ncm</p> <p>Scan post is torqued with max. 10 Ncm</p> |
| <p>Universal Casting Abutment Screws</p>  |  <p>T6 2024 Short (ratchet)</p> <p>T6 2025 Long (ratchet)</p> | <p>Max. 10 Ncm</p> |

|  T0 25900 | | |
|---|---|---------------|
| Products | Tighten Parts | Torque Values |
| Multi Unit Abutments  |  T6 2024 Short (ratchet) T6 2025 Long (ratchet) | 20 Ncm |
| Multi Unit Abutments (Straight/ 0°)  |  T0 25205 | Max. 30 Ncm |
| Ball Abutments  |  T0 25206 | Max. 30 Ncm |
| Equator® Abutments  |  774CHE Tightening Part (ratchet/screwdriver) | Max. 30 Ncm |

NucleOSS Dental Implant System - Implant Prospectus

Definitions

NucleOSS Dental Implants (T1, T3, Tpure, T5 Short, T6, T6 TorQ and Slimex) are medical devices that intended to be surgically placed in the bone of maxillary and/or mandibular arches to support prosthetic restorations for restoration of the patients' chewing function. The implants are manufactured out ISO 5832-2 pure titanium Grade 4 to fulfill the requirements of 93/42/EEC medical devices directive.

NucleOSS Dental Implant is recommended for delayed loading after 12 weeks. NucleOSS Implants are non-pyrogenik. Implants are provided as sterile in a hermetic tube with a blister package and are for single use. Do not sterilize again.

Surface Specification

All NucleOSS Dental Implants have **maxOcell®** surface. The surface is sand blasted and processed in acid baths. Following its special cleaning steps, the surface gains its hydrophilicity.

Patient Information

Surgery site must be clean before and after operation. The patient must be well informed about cleaning of the surgical area. The implant site should not be exposed with pressure due to chewing function. For a more detailed explanation please refer to Patient Information Brochure (D05.001/en).

Packaging

NucleOSS Dental Implants are protected with double barrier system as primary packaging (tube) and secondary packaging (blister). Implants are placed in blister package and are sterilized; and delivered to final packaging process as sterilized.

Labelling

NucleOSS Dental Implants are provided with tube label, tyvek label and box label. These labels enables to control the lot number in three different stage. The diameter label that is present on the tube states the diameter and length of the implant. Diameter labels are color coded for different diameters.

Sterilization

NucleOSS Dental Implants are provided as sterile. Implants are sterile washed under Class 10.000 and packed under Class 100 clean room technology. Implants are placed in hermetically sealed tubes and sterilized by 25-40 kGy of gamma irradiation. Gamma sterilization is protected until the best before if the packaging is not harmed. The proper storage conditions that is present on the package are 18-28° and %40-%60 moisture must taken seriously to protect sterilization.

⚠ Cautions

- NucleOSS Dental Implant System products are produced by Şanlilar Med. Kim. San. Tic. Ltd. Şti.. Brand usage and sale of the products can be conducted only by Şanlilar and partner firms.
- NucleOSS Dental Implant system models (T1, T3, Tpure, T5 Short, T6, T6 TorQ and Slimex) should apply accordingly with NucleOSS Surgical and Prosthesis Surgical Protocols. It is advised to wait at least 3 months (12 weeks) for osseointegration after surgery.
- Implants are in sterile packaging.
- Do not open the tube cover before use under any circumstances. Do not use the implant tubes that is opened or became deformed.
- Keep it away from the child reach.
- If the impression and transfer part is present inside the implant package, this part might broken in case of using it as tightening part. Tightening must be actualized with tightening part.
- It is important to decide on the implant that is suitable as diameter and length for implementation site. Advised implementation sites are present in the relevant product catalogs. Please refer to the relevant catalogs for further information (for T6 D01.001/en; for Tpure D01.002/en; for T3 D01.003/en; for T1 D01.004/en; for T6 TorQ D01.005/en; for Slimex D01.006/en and for T5 short D01.007/en document numbered.)
- Apply adequate number of implants with suitable diameters in an axis with compatible with the dentition.
- Use the proper drill compatible with diameter and length.
- Inform the patient before and after surgery.
- Adequate general state of health of the patient is necessary.
- Implant tube must removed from the blister to sterile laboratory cloth by wearing handgloves.

- Implants with broken protection ring or harmed package must not be used.
- To avoid the risk of contamination, the implant must insert to the bed in the moment of taking it out of the sterile tube without touching it or contacting anywhere.
- Bone development of the patient must be considered for the implant treatment.

Storage Conditions

Devices should be stored at 18°C - 28°C temperature range and 40% - 60%

Caution ! Previously used implants cannot be used again. Use only the NucleOSS Dental System components and surgical kits. Manufacturer will not take responsibility in case of using the parts other systems.

Caution ! If the tube consists of which is also colored accordingly with its platform diameters, must transfer on implant by hex driver and handtight to cover the implant. Otherwise, the tissue might fill the implant and create risk factor.

Traceability

A Lot number is written in each package. Also, each package contains three Lot number sticker. In order to trace back the product, this Lot number stickers must be attached to the patient's file and panoramic x-ray.

Indications

NucleOSS Dental Implants (T1, T3, Tpure, T5 Short, T6, T6 TorQ nad Slimex) are medical devices that intended to be surgically placed in the bone of maxillary and/or mandibular arches to support prosthetic restorations for restoration of the patient's chewing function.

Contraindications

Hypertension, cardiovascular diseases, diabetes bone metabolism disturbances, uncontrolled bleeding disorders, inadequate wound healing capacity, inadequate oral hygiene, serious internal medical problems, maxillary and mandibular growth not completed, poor general state of health, uncooperative and unmotivated patient, drug/alcohol/tobacco abuse, psychoses, long term treatment resistant functional disorders, xerostomia, granulocytopenia, Ehlers Danlos syndrome, radiation therapy, osteoradionecrosis, weak immune system, kidney failure, organ transplantation, fibrous dysplasia, crohn disease, use of steroids, corticosteroids or anticonvulsant usage, prophylactic antibiotics, creatinin, serum calsiun, titanium allergy, uncontrolled endocrine disorders, anticoagulation medicines, hemorrhagic diathesis, bruxism, parafunctional habits, unfavorable anatomic bone condition, uncontrolled periodontitis, temporomandibular joint disorders, treatable pathologic of jaw, changes in the oral mucosa, pregnancy, breast-feeding, osteoporotic crush fracture, respiratory disease, tyhroid or paratyhroid diseases. Active osteolytic patients, inflammatory, diagnosed malignancy, modular enlargements, tenderness, the patients with unexplained lump in the neck or head, infectious process in the implanting site, unrealistic patient expectations, unattainable prosthodontics reconstruction, lack of adequate training of pyhscian.

Risks

The risks may include inadvertent perforation of the nasal and maxillary sinus, local and systemic infections, perforation the soft tissue, nerve damage, temporary bumps and pain due to implantation, speech problems nad gingivitis.

Nerve, local or systemic bacterial infections and inactive endocarditis in susceptible individuals are included as the long term problems. Inaccurate implantation might risk the present dental axis.

NucleOSS Dental Implant System - Implant Prospectus

Surgical Guide

Do not insert the implant in high torque or speed (15 rpm). If it feels any difficulty, the practitioner must step back, extract the implant and check for the implant bed and drills.

MRI Safety Information

Non-clinical testings demonstrated that the NucleOSS Dental Implant System is MR Conditional. A patient with this device can be scanned safely, immediately after placement in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial field gradient of 4,000 -gauss/cm (40 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the NucleOSS Dental Implant System is expected to produce a maximum temperature rise of less than 1.7° after 15 minutes of continuous scanning. In non-clinical testing, the image artifact caused by the device extends approximately 5 mm from the implant when imaged with a gradient echo pulse sequence and a 3.0 MRI system.

Superstructures

NucleOSS Dental Implant System offers dentists a rich abutment options that can be used for all bodies with different purposes. These are defined as cemented solutions, screw-in solutions, custom solutions, overdenture solutions, cad-cam solutions and scanning bodies, auxiliary prosthetic parts, temporary parts. The product compability according to platform diameter and lenght of abutments can be seen in relevant catalogs of the products.

Abutments are only suitable for the bodies of NucleOSS Dental Implant System. The failure that is caused by using with other system might hurt the patient.

All Abutment belong to NucleOSS Dental Implant System are produced out ISO 5832-3 Grade-5 Titanium and ISO 5832-12 Cobalt-Chromium alloy and non-sterile and presented for sale. The must sterilize by the last user.

Sterilization

Titanium superstructures recommended sterilization condition as follows;

| Method | Moist heat sterilization according to ISO 17665 |
|---------------------|---|
| Cycle | Pre vacuum |
| Temperature | 132° C / 270° F |
| Exposure time | 4 minutes |
| Pressure | 2.2 Bar |
| Pre-vacuum | 3 times < 60 mbar |
| Minimum drying time | 30 minutes in chamber |

- Do not storage the products after sterilization.
- Minimum validated sterilization time and temprature required to achieve a -10⁻⁶ sterility assurance level.

Note: Sterilization parameters and methods shown are validated by NucleOSS Dental Implants. According to EN ISO 17665, the final responsibility for validation of sterilization rechniques and equipment lies directly with the practitioner. All autoclaves/sterilizers should be validated and maintained in accordance with EN ISO 17665-1.

Indications

NucleOSS abutments and prosthetic parts are intended for use with NucleOSS Dental Implant in the maxillar and mandibular arches to prived and support for crowns, bridges or overdentures for edentulous or partially edentulous cases.

Contraindications

Allergy may be developed to submaterials of the superstructures. Titanium, Cr-5 titanium alloy (titanium-aluminium-vanadium), PEEK (Polyetheretherketone), POM (Polyoxymethylen), ISO 5832-12 raw material CoCr alloy allergies must be considered.

Cautions

In case of not considering the cautions below, complications may occur as parts slipping into trachea or getting swallowed.

For dental technicians:

Make sure to protect the parts of abutments that goes inside the implant or this process must complete while abutment is connected to analog.

Abutment must be placed on the model by controlling if its in the right direction and must make sure that it fits perfectly. Abutment must be stabilized by abutment screw and hadtight only (max. 10 Ncm). After this point, the accurate prosthesis fir the treatment must be formed and extracted from the model by proper hex driver.

For dental professionals:

Clinician takes the abutments comes from the laboratory, and takes the temporary restorations, cover screws or gingiva formers out of the patient's mouth. After cleaning, disinfecting and sterilizing the prosthesis parts coming from the laboratory, the application must be followed as the given appropriate torque values.

- In case of tightening abutments more or less than the adviced torque values, abutment or implant might fail.
See below for recomendated torque values;

Torque Values for NucleOSS Parts

| | |
|---------------------------------------|-----------------------|
| All Covers and Gingiva Formers | Hand/Ratchet (10 Ncm) |
| All Cemented Abutment Screws | Max. 30 Ncm |
| Temporary (PEEK) Parts and Screws | 10 Ncm |
| Multi Unit Abutments | Max. 30 Ncm |
| Abutment Occlusional Screw | 20 Ncm |
| All Cap Screws | 10 Ncm |
| All Universal Casting Abutment Screws | Max. 30 Ncm |
| Ball Abutments | Max. 30 Ncm |
| Equator® Abutments | Max. 30 Ncm |
| CAD-CAM Ti-Base Abutment Screws | Max. 30 Ncm |

Caution: Products must be used right after sterilization. Sterilization conditions are given above.

Covering the Screw Channel:

Before inserting the crown on abutment, abutment screw channel must be covered with sealing component (teflon). This process enables to take the abutment out if that is nedeed afterwards.

NucleOSS Dental Implant System - Implant Prospectus

Adjustment to the Patient's Anatomy:

PEEK, Titanium or Titanium Alloy abutments may shortened until the top level of the connection screw, if needed to adjust to the patient's anatomy.

CoCr Castable Abutment Cautions:

CoCr raw material must comply with ISO 5832-12 standard. The molten metal temperature should not exceed 1420 °C to avoid the melting of CoCr Casting Abutment.

CoCr Coefficient of Thermal Expansion (CTE) (α):

- For 20 - 100°C, 13.2 (10-6 K-1)
- For 20 - 200°C, 13.3 (10-6 K-1)
- For 20 - 300°C, 13.5 (10-6 K-1)
- For 20 - 400°C, 13.8 (10-6 K-1).

Gingiva height should not exceed 5.0 mm and abutment angle should not exceed 25°.






















Caution! Temporary abutments should not stay in patients' mouth more than exposure time. Exposure time for temporary abutment (Ti, PEEK) is 180 days.

-Şanlılar Tıbbi Cihazlar Med. Kim. San. Tic. Ltd. Şti. consistently gives trainings. While preparing this prospectus, it is assumed that the clinician participates in these trainings and have prior knowledge to the practise. For further information please refer to surgical manuals (for T6 D03.001/en, for Tpure D03.002/en, for T3 D03.003/en, for T1 D03.004/en, for T6 TorQ D03.005/en, for Slimex D03.006/en and for T5 Short D03.007/en document numbered) and surgical instruments instructions for use (D02.002 document numbered).

Any serious incident related to the medical device should be reported to the manufacturer and the local authority.

If the storage conditions cannot be met, the product should not be used.

Warning Signs

| | | | | | |
|--|--|---|--------------------------------|--|--------------------------------|
|  REF | Catalog number |  i | Consult operating instructions |  MD | Medical Device |
|  LOT | Lot number |  Hourglass | Use before expiry date |  CC | Country of Manufacturer |
|  Single use | For single use only |  Factory | Manufacturer |  UDI | Unique Device Identifier |
|  STERILE | Sterile by gamma irradiation (for implant) |  Rx only | For prescription use only |  # | Model number |
|  Do not use if package is damaged | Do not use if package is damaged |  MR | MRI Safety Information |  Do not sterilize (for implant) | Do not sterilize (for implant) |
|  Temperature limit | Temperature limit |  Cautions | Cautions | | |
|  Humidity limit | Humidity limit |  Non Sterile (For Abutment) | Non Sterile (For Abutment) | | |
|  EC REP | Authorised Representative |  Double Sterile Barrier System | Double Sterile Barrier System | | |



Şanlılar Tıbbi Cihazlar Medikal Kimya San. Tic. Ltd. Şti.
İTOB Organize Sanayi Bölgesi 9. Sk. No:24 D:26, 35865 Menderes / İzmir / Türkiye
T: 90 232 799 03 04 (pbx) F: +90 232 799 03 06

www.NucleOSS.com



EC REP Nucleoss Europe GmbH
Graben 17, 64646 Heppenheim (Bergstr.)
GERMANY

Notes



Notes



Notes



Notes



Şanlılar Tıbbi Cihazlar
Medikal Kimya San.
Tic. Ltd. Şti



10018 Sk. No:7 ITOB Organize
Sanayi Bölgesi 35477,
Tekeli - Menderes - İZMİR

T +90 232 799 03 04
F +90 232 799 03 06
E info@nucleoss.com



NucleOSS Europe GmbH
Craben 17
64646 Heppenheim / Germany

T +49 (0) 6252 795 7772
F +49 (0) 6252 795 7773
E europe@nucleoss.com

innovative **vision**

