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Assoc. Prof. Dr. Ismail Doruk Kocyigit graduated from Ankara University Dentistry Faculty in 2002 and in the same year he started his residency training in Oral and Maxillofacial Surgery at the same faculty. He has been assigned Assoc. Prof. Dr. at Kirikkale University, Department of Oral and Maxillofacial Surgery in 2014 and he is currently at the same position in the same university. He is a member of (IAOMS), (EACMS), (ACBID), (CAI Academy). He was joined an observer ship programme in London 2012 with the support of ACBID Scholarship. His principle areas of clinical practice and research are and Reconstruction, Orthognatic Surgery, Cleft Lip and Palate, Trauma&Dental Implants. Dr. Kocyigit has 32 published/ahead of publishing papers in SCI/SCI-expanded journals and severals in national journals and 5 International Scientific Awards. He has been several national and international lectures about his interested topics.

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Patient Information:


46 years old non-smoker male patient with no systemic disease.

Patient History: The patient referred to clinic for the treatment of missing 16 and 17.

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The fixed restoration treatment was decided for the treatment of the missing 16 and 17 teeth.

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The patient was informed that sinus augmentation procedure needed to place implant due to resorption, inadequate bone around maxillary sinus cavity. After the patient's accepting the treatment, external sinus lift operation be carried out.

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Allogenic graft material and resorbable collagen membrane are used during external sinus lifting operation. The healing period was 8 months. Post-operative radiography showed adequate bone healing on the operated site.

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The site is re-opened after 8 months and osteotomized for Nucleoss T6 implants. On the left, osteotomy and healed sinus wall can be seen.

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For missing tooth number 16 4,1 mm diameter and 10 mm of length; for missing tooth 17 4,8 mm diametered and again 10 mm of length Nucleoss T6 implants are placed.

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The radiography that was taken right after the operation is on the left.

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The intraoral view of the cemented suprastructure after 44 months clinical follow-up.

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Periapical radiography after 44 months clinical follow up.