M.Sc.Dr. Khagani Huseynov

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He graduated from Azerbaijan Medical University Faculty of Dentistry in 2001. In 2006, he completed the expertise course on Maxillofacial Surgery and Implantology in Post-Diploma Training Academy, St. Petersburg University. In 2011, he received his master degree in laser dentistry at Achen University, Germany. In 2012, he completed the surgical courses of Prof. Dr. Ateş Parlar. In 2015 he graduated from the education program of Implantology in Frankfurt am Main University, Germany. He is also the head of the Azerbaijan Esthetics Dentistry Academic Association). He has been various organising congresses and giving seminars and courses in his country since 2008.

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
Patient Information

50 years old, with poor mouth hygiene, male. The patient is smoking 20 cigarettes per day.

Patient History:

The patient consulted to clinic with the complaint of changing color in the front tooth and loose tooth.

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
Due to trauma accrued many years ago, the root of the tooth 21 was subject to resorption. In addition to this, the tooth itself grew, rotated, foregrounded and placed itself in the direction of mesial.

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The tooth has extracted and implant (Nucleoss T4 D5.0mm/L 14mm) placed in the direction of palatinal.

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The implant was covered with xenograft bone material. Gingiva former was placed to implant with collagen sponge(membran).

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After a day of taking impression and preparing the temporary prosthetics, its applied in a way that there will not be a functional load.

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The visual was taken 4 months ago. The situation of the soft tissue composed around the implant and emergence profil.

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Prepared abutment attached to implant with 25N/cm torque.

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The visual of the metal-ceramic in mouth. For esthetic reasons, fully resemblance with the rest of the teeth provided.

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Mounting the implant in the exact position and by means of bone addition protecting the vestibular wall fully.

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Post operation x-ray visual. Implementation places and follow-up periods; In the teeth number 46 - T3 Ø4.2mm/L12mm (7 years), in the tooth 27 - T3 Ø5.0mm/L14mm (7 years), in the teeth 32 and 42 -T2 Ø3.8mm/L14mm (8 years 7 months), in the tooth number 21 - T4 Ø5.0mm/L14mm (2 years)

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Intra-oral image of the case; 9 years after implementation of T2 and 2 years after implementation of T4 implants.